

## **EXHIBIT 16**

(7)

Consecutive digestion of APTag/BglII + HindII

BglII digestion 1.9 ml

water 1.7 ml

React I 2 ml

HindII 2 ml

40 ml at 37°C

(purify by elutip method)

05:00:06

digestion of 4-18B PCR products

4-18B (#3, #4) 16 ml 4-18B (5, 5') 16 ml

React 3 2 ml React 3 2 ml

spermidine 1 ml spermidine 1 ml

EcoRI 1 ml BglII 1 ml

20 ml 20 ml

~~Phenol/chloroform ext. (no test)~~ at 37°C 05:00:06


add React 1 2 ml React 1 2 ml

water 16 ml water 16 ml

StuI 2 ml HindII 2 ml

40 ml 40 ml

06:45



## Ligation

	Exp.	Control (self)
pPR998/ <del>XmnI</del> RI (0.5 μg/μl)	4.0 μl	4.0 μl
4-180 PCR/ <del>SmaI</del> RI (12 μg/μl)	0.5 μl	—
5X buffer	2.0 μl	2.0 μl
water	3.0	3.5
ligase	0.5 μl	0.5
	10.0 μl	10.0

Legation

APtag / Hind III - B3/II (150 ng /  $\mu$ l)4-1 BB (PCR) / Hind III - B3/II (80 ng /  $\mu$ l)

5X buffer

water

ligase

exp.

1  $\mu$ l2  $\mu$ l2  $\mu$ l

4.5

0.5

10  $\mu$ l

self

1  $\mu$ l

-

2  $\mu$ l

6.5

0.5  $\mu$ l120  $\mu$ l

method (self)

4.0  $\mu$ l

-

2.0  $\mu$ l

3.5

0.5

10.0

